# medivon.

User manual

PL | EN | DE | ES | IT | FR | CZ  
SL | PT | LT | RU | SK

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Start by selecting the appropriate mode: heating or cooling. Make sure the head is securely attached and not damaged. Then turn on the device by holding down the button corresponding to the mode visible on the bottom of the head for a while. After starting, pressing the button cyclically changes the heating or cooling power on a 2-step scale. To turn off the heads, hold down the mode button again for a few seconds.

Note: we replace the massage heads only when the device is turned off. It is enough to remove them from the device by pulling (they do not have a screw thread).

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The manufacturer declares under his sole responsibility that the Medivon Reva device complies with Directive 2014/30/EU of the European Parliament and of the Council of February 26, 2014 on the harmonization of the laws of the Member States relating to electromagnetic compatibility, as well as Directive 2011 of the European Parliament and of the Council /65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to the Directive of the European Parliament and 2011/65/EU as regards the list of restricted substances. The full text of the declaration of conformity is available at the internet address

Importer: Armare S.A. | ul Grzybowska 87, Warszawa, 00-844, Polska  
Model: Medivon Reva

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Electrical and electronic devices as well as batteries and accumulators must not be mixed with other household waste. This is indicated by the symbol of the crossed out bin placed on the packaging. Worn-out electric and electronic equipment as well as batteries and accumulators should be delivered to local separate waste collection points or to the seller, and for details inquire in your commune. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution during further disposal or recycling of the equipment.